

SonoTouch 30







PROFESSIONAL Clinical Application UI


The dedicated icon-based presets for major MSK / Anaesthesia/ Vascular positions will help doctors acquire optimal images in an easier, quicker, and more intuitive way.

It uses the latest Supper Tissue technology on the 11MHz Linear probe to enhance the near field image within 2.5 cm depth which is the focal region of most Point of care applications.











Patella Tendon Tibian Insertion Part, B Mode




Palmaris Longus Tendon, B Mode





Liver ,B Mode




Carotid, C Mode




Central Venous Access, B Mode





Popliteal Sciatic Nerve, B Mode



Ulnar Artery ,PW Mode



Kidney ,B Mode

Super Needle (Option)

Super Needle(Option) is our needle enhancement technology supporting angle adjustment up to +/- 20 degrees. It enhances the needle's appearance for a more confident approach to anatomical targets.



Without Super Needle Needle is Not clear



With Super Needle Needle is Better Visualized











WWW.CHISON.COM Trolley(TR-1) Tri-connectors

Standard Configuration

Main unit,10.4"Touch screen,1 probe connector, 3 USB ports,VGA port,LAN port

Options



5MHz-10MHz Linear



5MHz-11MHz Linear



2.5MHz-5MHz Convex C3






2.5MHz-4MHz Micro-Convex 2.5MHz-4MHz Phased Array MC3 Р3



4MHz-7.5MHz Transvaginal



5MHz-10MHz Trans-Rectal



4MHz-8MHz Micro-Convex



4.5MHz-7MHz Micro-Convex MC5V

- Biopsy kit is available (probe dependent)
- Video printer, PC printer
- Super Needle
- Trolley(TR-1)
- DICOM3.0
- Tri-connectors
- Suitcase(CS-1)



FULL SCREEN display mode

(Height adjustable)

CHISON Medical Technologies Co., Ltd.

Sales & Service Contact Address:

No.9, Xinhuihuan Road, Xinwu District, Wuxi, Jiangsu, China 214028

TEL: 0086-510-85310593 / 85310937

FAX: 0086-510-85310726 EMAIL: export @chison.com.cn